

**Livsmedel – Bestämning av saxitoxin och
dc-saxitoxin i musslor – HPLC-metod med
postkolonn derivatisering**

**Foodstuffs – Determination of saxitoxin and
dc-saxitoxin in mussels – HPLC method using
post column derivatisation**

Den europeiska förstandarden ENV 14194:2002 gäller som svensk standard. Detta dokument innehåller den officiella engelska versionen av ENV 14194:2002.

The European Prestandard ENV 14194:2002 has the status of a Swedish Standard. This document contains the official English version of ENV 14194:2002.

Dokumentet består av 13 sidor.

Upplysningar om **sakinnehållet** i standarden lämnas av SIS, Swedish Standards Institute, tel 08 - 555 520 00.

Standarder kan beställas hos SIS Förlag AB som även lämnar **allmänna upplysningar** om svensk och utländsk standard.

Postadress: SIS Förlag AB, 118 80 STOCKHOLM
Telefon: 08 - 555 523 10. *Telefax:* 08 - 555 523 11
E-post: sis.sales@sis.se. *Internet:* www.sisforlag.se

EUROPEAN PRESTANDARD
PRÉNORME EUROPÉENNE
EUROPÄISCHE VORNORM

ENV 14194

April 2002

ICS 67.120.30

English version

Foodstuffs - Determination of saxitoxin and dc-saxitoxin in mussels - HPLC method using post column derivatisation

Produits alimentaires - Détermination de la teneur en saxitoxine et en dc-saxitoxine dans les moules - Méthode par chromatographie liquide haute performance après dérivation post-colonne

Lebensmittel - Bestimmung von Saxitoxin und DC-Saxitoxin in Muscheln - HPLC-Verfahren mit Nachsäulenderivatisierung

This European Prestandard (ENV) was approved by CEN on 14 December 2001 as a prospective standard for provisional application.

The period of validity of this ENV is limited initially to three years. After two years the members of CEN will be requested to submit their comments, particularly on the question whether the ENV can be converted into a European Standard.

CEN members are required to announce the existence of this ENV in the same way as for an EN and to make the ENV available promptly at national level in an appropriate form. It is permissible to keep conflicting national standards in force (in parallel to the ENV) until the final decision about the possible conversion of the ENV into an EN is reached.

CEN members are the national standards bodies of Austria, Belgium, Czech Republic, Denmark, Finland, France, Germany, Greece, Iceland, Ireland, Italy, Luxembourg, Malta, Netherlands, Norway, Portugal, Spain, Sweden, Switzerland and United Kingdom.



EUROPEAN COMMITTEE FOR STANDARDIZATION
COMITÉ EUROPÉEN DE NORMALISATION
EUROPÄISCHES KOMITEE FÜR NORMUNG

Management Centre: rue de Stassart, 36 B-1050 Brussels

Foreword

This document (ENV 14194:2002) has been prepared by Technical Committee CEN /TC 275, "Food analysis - Horizontal methods", the secretariat of which is held by DIN.

According to the CEN/CENELEC Internal Regulations, the national standards organizations of the following countries are bound to announce this European Prestandard: Austria, Belgium, Czech Republic, Denmark, Finland, France, Germany, Greece, Iceland, Ireland, Italy, Luxembourg, Malta, Netherlands, Norway, Portugal, Spain, Sweden, Switzerland and the United Kingdom.

1 Scope

This European Prestandard specifies a method for the quantitative determination of saxitoxin (STX), dc-saxitoxin (dc-STX) and the qualitative determination of neo-saxitoxin, and the gonyautoxins GTX-2 and GTX-3 in mussels. The method can also be used to identify the N-sulfocarbamoyl toxins C-1, C-2, GTX-5 and GTX-6 after hydrolysis and, if these toxins are present, to exclude false positive results for GTX-2, GTX-3, neo-saxitoxin and saxitoxin. For mussel the limit of quantification is for saxitoxin 0,04 mg/kg mussel meat and for dc-saxitoxin 0,03 mg/kg mussel meat (signal/noise = 10).

The limits of detection for C-1, C-2, GTX-2, GTX-3, GTX-5, GTX-6 and neo-saxitoxin have not been determined.

2 Normative references

This European Prestandard incorporates by dated or undated reference, provisions from other publications. These normative references are cited at the appropriate places in the text and the publications are listed hereafter. For dated references, subsequent amendments to or revisions of any of these publications apply to this European Prestandard only when incorporated in it by amendment or revision. For undated references the latest edition of the publication referred to applies (including amendments).

EN ISO 3696:1995, *Water for analytical laboratory use - Specification and test methods (ISO 3696:1987)*

3 Principle

PSP (Paralytic Shellfish Poisoning) toxins are extracted from mussels homogenate with an acidic aqueous solution. After centrifugation the supernatant is purified by solid phase extraction (SPE) over a C18 clean-up cartridge. Part of the extract is directly injected on the HPLC for toxin determination. After the separation on the analytical column post column derivatisation is used, with oxidation with periodic acid and after acidification the derivatized toxins are detected using fluorimetric detection. HPLC analysis is separated in two parts: a system to separate the saxitoxin group toxins (saxitoxin, dc-saxitoxin and neo-saxitoxin) and a system to separate the GTX group toxins (GTX-2 and GTX-3). Identification after hydrolysis is based on the transformation of the N-sulfo-carbamoyl toxins to their corresponding carbamoyl toxins.

WARNING - PSP toxins are strong neurotoxins. Gloves and safety glasses should be worn at all times, and all standard and sample preparation stages should be carried out in a fume cupboard.

4 Reagents

4.1 General

During the analysis, unless otherwise stated, use only water according to grade 1 of EN ISO 3696:1995.

All chemicals shall be of pro analysis (p.a) quality, unless otherwise indicated.

Reference materials (calibrants of the toxins) originating from other sources than indicated may also be used if well-characterised and with a well-defined mass concentration.

4.2 Methanol

4.3 Acetonitrile, HPLC quality

4.4 Hydrochloric acid, volume concentration $\varphi(\text{HCl}) \approx 25 \%$, (acidimetric)

4.4.1 Hydrochloric acid solution, substance concentration $c = 1,0 \text{ mol/l}$

Dilute 130 ml of hydrochloric acid solution (4.4) with 870 ml of water.

4.4.2 Hydrochloric acid solution, $c = 0,1 \text{ mol/l}$

Dilute 100 ml of hydrochloric acid solution (4.4.1) with 900 ml of water.

4.5 Octane sulfonic acid solution, $c = 0,5 \text{ mol/l}$

Dissolve 2,9 g of octane sulfonic acid sodium salt monohydrate in 25,0 ml of water.

4.6 Phosphoric acid, the mass fraction is 85 %

4.6.1 Phosphoric acid solution, $c = 0,5 \text{ mol/l}$

Take 6,7 ml of phosphoric acid (4.6) and dilute to 200,0 ml with water.

4.7 Ammonia solution, $\varphi = 25 \%$ and $\varphi = 1 \%$

For the 1 % solution, dilute 40 ml ammonia solution $\varphi = 25 \%$ with 960 ml of water.

4.8 Periodic acid solution, $c = 0,01 \text{ mol/l}$

Dissolve 1,14 g of periodic acid in water and dilute up to 500 ml with water.

4.9 Acetic acid, 100 %

4.9.1 Acetic acid solution 1, $c = 1 \text{ mol/l}$

Take 57,2 ml of acetic acid solution (4.9) and dilute to 1,0 l with water.

4.9.2 Acetic acid solution 2, $c = 0,2 \text{ mol/l}$

Dilute 200 ml of acetic acid solution 1 (4.9.1) with 800 ml of water.

4.9.3 Acetic acid solution 3, $c = 0,03 \text{ mol/l}$

Dilute 150 ml of acetic acid solution 2 (4.9.2) with 850 ml of water.

ENV 14194:2002 (E)

4.10 Helium**4.11 Saxitoxin standard solutions****4.11.1 Saxitoxin stock solution, mass concentration $\rho = 1,0 \mu\text{g/ml}$**

Prepare a calibration solution containing $1,0 \mu\text{g/ml}$ of saxitoxin in acetic acid solution 3 (4.9.3). Store the calibration solution cool (at approximately $+ 4 \text{ }^\circ\text{C}$) and in the dark. This solution is stable for at least 6 months.

NOTE 1 An ampoule containing approximately $0,2 \text{ ml}$ of saxitoxin calibration solution, mass concentration of saxitoxin of $0,14 \text{ mg/ml}$ is for example available from the National Research Council Canada, Halifax, Canada (set of calibrants PSP-1). In order to obtain a calibration solution of $1,0 \mu\text{g/ml}$, the content of the ampoule can be quantitatively transferred to $27,8 \text{ ml}$ acetic acid 3 (4.9.3) (total volume $28,0 \text{ ml}$). The total volume can be checked and adjusted by weighing (5.9), using density of water for calculations. For this purpose weigh the ampoule before opening, take care that all glass is preserved after opening and weigh again after the calibrant solution is transferred.

NOTE 2 A calibration solution containing GTX-1 and GTX-4 is also available and may be used for qualitative determination of these toxins. This has not been tested in the certification study of the Standard Measurements & Testing Programme, designed to certify the mass fractions of saxitoxin and dc-saxitoxin in mussel reference materials, see [6].

4.11.2 Saxitoxin calibration solutions, $\rho = 0,05 \mu\text{g/ml}$ to $0,30 \mu\text{g/ml}$

Prepare for the determination of saxitoxin a calibration series of increasing concentration in the mass concentration range of $0,05 \mu\text{g/ml}$ to $0,30 \mu\text{g/ml}$ from saxitoxin stock solution (4.11.1) in acetic acid solution 3 (4.9.3). Prepare fresh every day of analysis.

4.12 Dc-saxitoxin standard substance

Store cool and in the dark.

4.13 Dc-saxitoxin stock solution, $\rho = 1,0 \mu\text{g/ml}$

Prepare from dc-saxitoxin standard substance (4.12) a dc-saxitoxin calibration solution $1,0 \mu\text{g/ml}$ in acetic acid solution 3 (4.9.3) by weighing (5.9), using density of water for calculations. This solution is stable for at least 5 months (see Bibliography, [3]).

4.14 Dc-saxitoxin calibration solution, $\rho = 0,05 \mu\text{g/ml}$ to $0,30 \mu\text{g/ml}$

Prepare for the determination of dc-saxitoxin a calibration series of increasing concentration in the mass concentration range of $0,05 \mu\text{g/ml}$ to $0,30 \mu\text{g/ml}$ from Dc-saxitoxin stock solution in acetic acid solution 3 (4.9.3). Prepare fresh every day of analysis.

NOTE For the quantitative determination of saxitoxin and dc-saxitoxin a combined calibration solution for these 2 toxins can be made if the valley between the 2 toxin peaks in the HPLC chromatogram is $\leq 10\%$ of the sum of the peak heights. If this condition is fulfilled in the HPLC separation, one combined calibration solution containing saxitoxin and dc-saxitoxin can be prepared in the mass concentration range of $0,05 \mu\text{g/ml}$ to $0,30 \mu\text{g/ml}$ (for both toxins) in acetic acid solution 3 (4.9.3).

4.15 Neo-Saxitoxin standard solutions**4.15.1 Neo-saxitoxin calibration solution, $\rho = 2,0 \mu\text{g/ml}$**

Prepare a calibration solution containing $2,0 \mu\text{g/ml}$ of neo-saxitoxin in acetic acid solution 3 (4.9.3). Store the calibration solution cool (at approximately $+ 4 \text{ }^\circ\text{C}$) and in the dark. This solution is stable for at least 6 months.

NOTE 1 An ampoule containing $0,2 \text{ ml}$ neo-saxitoxin calibration solution, mass concentration $0,14 \text{ mg/ml}$ is for example available from the National Research Council Canada, Halifax, Canada (set of calibrants PSP-1). In order to obtain a calibration solution of $2,0 \mu\text{g/ml}$, the content of the ampoule can be quantitatively transferred to $13,8 \text{ ml}$ acetic acid solution 3 (4.9.3) (total volume $14,0 \text{ ml}$). The total volume can be checked and adjusted by weighing (5.9), using density of water for calculations.