

Chemical disinfectants – Quantitative suspension test for the evaluation of virucidal activity against bacteriophages of chemical disinfectants used in food and industrial areas – Test method and requirements (phase 2, step 1)

ICS 11.080.20; 67.050; 71.100.35

Språk: engelska

Publicerad: augusti 2003

Europastandarden EN 13610:2002 gäller som svensk standard. Detta dokument innehåller den officiella engelska versionen av EN 13610:2002.

The European Standard EN 13610:2002 has the status of a Swedish Standard. This document contains the official English version of EN 13610:2002.

Dokumentet består av 46 sidor.

Upplysningar om **sakinnehållet** i standarden lämnas av SIS, Swedish Standards Institute, tel 08 - 555 520 00.

Standarder kan beställas hos SIS Förlag AB som även lämnar **allmänna upplysningar** om svensk och utländsk standard.

Postadress: SIS Förlag AB, 118 80 STOCKHOLM
Telefon: 08 - 555 523 10. *Telefax:* 08 - 555 523 11
E-post: sis.sales@sis.se. *Internet:* www.sis.se

EUROPEAN STANDARD
NORME EUROPÉENNE
EUROPÄISCHE NORM

EN 13610

December 2002

ICS 11.080.20; 67.050; 71.100.35

English version

Chemical disinfectants - Quantitative suspension test for the evaluation of virucidal activity against bacteriophages of chemical disinfectants used in food and industrial areas - Test method and requirements (phase 2, step 1)

Désinfectants chimiques - Essai quantitatif de suspension pour l'évaluation de l'activité virucide contre les bactériophages des désinfectants chimiques utilisés dans le domaine de l'agro-alimentaire et dans l'industrie - Méthode d'essai et exigences (phase 2, étape 1)

Chemische Desinfektionsmittel - Quantitativer Suspensionsversuch zur Bestimmung der viruziden Wirkung gegenüber Bakteriophagen von chemischen Desinfektionsmitteln in den Bereichen Lebensmittel, und Industrie - Prüfverfahren und Anforderung (Phase 2, Stufe 1)

This European Standard was approved by CEN on 23 October 2002.

CEN members are bound to comply with the CEN/CENELEC Internal Regulations which stipulate the conditions for giving this European Standard the status of a national standard without any alteration. Up-to-date lists and bibliographical references concerning such national standards may be obtained on application to the Management Centre or to any CEN member.

This European Standard exists in three official versions (English, French, German). A version in any other language made by translation under the responsibility of a CEN member into its own language and notified to the Management Centre has the same status as the official versions.

CEN members are the national standards bodies of Austria, Belgium, Czech Republic, Denmark, Finland, France, Germany, Greece, Iceland, Ireland, Italy, Luxembourg, Malta, Netherlands, Norway, Portugal, Spain, Sweden, Switzerland and United Kingdom.



EUROPEAN COMMITTEE FOR STANDARDIZATION
COMITÉ EUROPÉEN DE NORMALISATION
EUROPÄISCHES KOMITEE FÜR NORMUNG

Management Centre: rue de Stassart, 36 B-1050 Brussels

Contents

	page
Foreword.....	3
Introduction	4
1 Scope	5
2 Normative references	6
3 Terms and definitions.....	6
4 Requirements	7
5 Test methods.....	7
5.1 Principle	7
5.2 Material and reagents	7
5.3 Apparatus and glassware	12
5.4 Preparation of host bacteria suspensions	13
5.5 Preparation of bacteriophage suspension.....	14
5.6 Product test solution	16
5.7 Procedure	16
5.8 Calculation and expression of results.....	19
5.9 Test report	23
Annex A (normative) Test for validation of dilution-neutralization and molecular sieving methods.....	25
Annex B (informative) Neutralizers.....	33
Annex C (informative) Example of a typical test report.....	35
Annex D (informative) Information on the application and interpretation of European standards on chemical disinfectants and antiseptics	39
Annex E (informative) Example of plaques from lysates of phages P001 and P008.....	41
Bibliography	44

Foreword

This document (EN 13610:2002) has been prepared by Technical Committee CEN /TC 216, "Antiseptics and chemical disinfectants" the secretariat of which is held by AFNOR.

This European Standard shall be given the status of a national standard, either by publication of an identical text or by endorsement, at the latest by June 2003, and conflicting national standards shall be withdrawn at the latest by June 2003.

In this European Standard the Annex A is normative and the Annexes B, C, D and E are informative.

According to the CEN/CENELEC Internal Regulations, the national standards organizations of the following countries are bound to implement this European Standard: Austria, Belgium, Czech Republic, Denmark, Finland, France, Germany, Greece, Iceland, Ireland, Italy, Luxembourg, Malta, Netherlands, Norway, Portugal, Spain, Sweden, Switzerland and the United Kingdom.

Introduction

This European Standard describes a suspension test method intended to establish whether a product proposed as a disinfectant in the fields described in clause 1 has or does not have virucidal activity against bacteriophages.

NOTE Virulent bacteriophages (phages) lytic for (i.e. virulent for) starter cultures of *Lactococcus lactis* subsp. *lactis* used for the production of cheese and other fermented milk products are used as model viruses to test the virucidal activity of a product.

The laboratory test closely simulates practical conditions of application. Chosen conditions (contact time, temperature, viruses [i.e. bacteriophages] in suspension, ...) reflect parameters which are found in practical situations including conditions which may influence the action of disinfectants.

Each utilization concentration found from this test corresponds to defined experimental conditions.

The conditions are intended to cover general purposes and to allow reference between laboratories and product types.

However for some applications the recommendations of use may differ and therefore additional test conditions need to be used.

1 Scope

This European Standard specifies a test method (phase 2, step 1) and requirements for the minimum virucidal activity against bacteriophages of chemical disinfectants that form a homogeneous, physically stable preparation in hard water and that are used in food and industrial areas, excluding areas and situations where disinfection is medically indicated and excluding products used on living tissues.

This European Standard applies at least to the following:

a) processing, distribution and retailing of:

1) food of animal origin:

- milk and milk products;
- meat and meat products;
- fish, seafood, and their products;
- eggs and egg products;
- animal feeds;
- etc.;

2) food of vegetable origin:

- beverages;
- fruits, vegetables and their derivatives (including sugar, distillery ...);
- flour, milling and baking;
- animal feeds;
-;

b) other industrial areas:

- biotechnology (yeast, proteins, enzymes, ...).

Using this European Standard, it is not possible to determine the virucidal activity against bacteriophages of undiluted product as some dilution is always produced by adding the inoculum and interfering substance.

For chemical disinfectants that can be used without dilution it is not possible to determine whether these products, at a concentration above 80 % have a virucidal activity against bacteriophages.

NOTE The method described is intended to determine the activity of commercial formulations or active substances on viruses (bacteriophages) in the conditions in which they are used.

This European Standard is applicable only to disinfectants, complying with the validation test (see annex A).

EN 13610:2002 (E)**2 Normative references**

Not applicable.

3 Terms and definitions

For the purposes of this European Standard, the following terms and definitions apply.

3.1**product (for chemical disinfection and/or antiseptis)**

chemical agent or formulation used as a chemical disinfectant or antiseptic

[EN 1040:1997]

3.2**virucide against bacteriophages**

product which inactivates the bacteriophages

NOTE The adjective derived from "virucide" is "virucidal".

3.3**virucidal activity against bacteriophages**

capability of a product to reduce the infectivity of intact bacteriophages belonging to reference bacteriophage strains P001 and P008 for at least 4 lg under the conditions defined by this standard

3.4**infectivity of bacteriophages**

the ability of a bacteriophage to propagate in a suitable host bacterial cell resulting in the release of bacteriophage progeny

3.5**inactivation of bacteriophage**

the reduction of infectivity of a bacteriophage by a product specified as a chemical disinfectant

3.6**reference bacteriophage suspension**

bacteriophage suspension of a defined virus strain maintained in reference centers and which should not be passaged more than 10 times

3.7**stock bacteriophage suspension**

bacteriophage suspension of a defined strain that has been multiplied on a large scale to obtain a bacteriophage suspension revealing identical characteristics as the reference bacteriophage suspension

NOTE This stock bacteriophage suspension is used to prepare a high-titer test bacteriophage suspension for the bacteriophage inactivation test.

3.8**high-titer bacteriophage suspension**

the high-titer bacteriophages suspension obtained from agar plates revealing confluent lysis in the bacterial lawn that is used to prepare the bacteriophage test suspension in the virucidal testing of the disinfectant

3.9**bacteriophage test suspension**

the bacteriophage suspension of a defined titer that is used in the virucidal testing of the disinfectant

4 Requirements

The product diluted in hard water when tested in accordance with clause 5 shall demonstrate at least a 4 log₁₀ reduction of infectivity of bacteriophages when tested in the presence of a volume fraction of 1 % acidic whey (prepared from acidified low-fat milk) or optionally in the presence of a volume fraction of 1 % skim milk as the interfering substance according to its practical applications and under the required test conditions (20 °C, 15 min, 2 reference bacteriophage strains).

The virucidal activity against bacteriophages shall be evaluated using the two virulent bacteriophages *Lactococcus lactis* subsp. *lactis* bacteriophage P001 and *Lactococcus lactis* subsp. *lactis* bacteriophage P008.

Both phages shall be propagated on the host strain *Lactococcus lactis* subsp. *lactis* F7/2.

The determined virucidal concentration of the tested product is suggested as being suitable for practical situations of use.

Where appropriate, additional specific virucidal activity against bacteriophages shall be determined under other conditions of time, temperature, additional strains and interfering substances in accordance with (see 5.7.1) in order to take into account intended specific use conditions.

NOTE For these additional conditions, the concentration defined as a result can be lower than the one obtained under the initial test conditions of 20 °C, 15 min, 2 selected bacteriophage reference strains.

5 Test methods

5.1 Principle

5.1.1 A test suspension of bacteriophages in a solution of interfering substances is added to a prepared sample of the product under test diluted in hard water.

The mixture is maintained at 20 °C ± 1 °C for 15 min ± 10 s (required obligatory test conditions).

At the end of the contact time, aliquots are taken and the virucidal activity against bacteriophages in this portion is immediately neutralized or suppressed by a validated method. The method of choice is dilution-neutralization with a validated neutralizer. If a suitable neutralizer is not available for a specific product, removal of the product by molecular sieving (i.e. gel filtration) shall be used.

The number of surviving bacteriophage particles and the number of bacteriophage particles in the test suspension are determined from appropriate dilution series with a factor of 10 prepared in medium in test tubes.

5.1.2 Additional and optional exposure times, temperatures and interfering substances are specified (see 5.7.1).

NOTE 1 The test described is based on an assessment (under specific conditions) which gives a reduction of at least 99,99 % (4 lg) of the infectivity of the different phages after different contact times.

NOTE 2 For principal reasons, the result of an inactivation applied to a viral population is not necessarily equal to 100 %: i.e., one cannot conclude that there is a 100 % inactivation when on conducting the experiment no infectious phage are found within a limit number of sampling.

5.2 Material and reagents

5.2.1 Test organisms

The virucidal activity against bacteriophages shall be evaluated using the two following bacteriophage strains :

EN 13610:2002 (E)

— *Lactococcus lactis* subsp. *lactis* bacteriophage P001 DSM 4262¹⁾;

— *Lactococcus lactis* subsp. *lactis* bacteriophage P008 DSM 10567.

Both phages shall be propagated on the following host strain:

— *Lactococcus lactis* subsp. *lactis* F7/2 DSM 4366.

5.2.2 Culture media and reagents

5.2.2.1 General

The reagents shall be of analytical grade and/or appropriate for microbiological purposes. They shall be free from substances that cause toxic or inactivating effects either to the bacteriophages or to the host bacteria used for counting phage-derived plaques.

NOTE To improve reproducibility, it is recommended that commercially available dehydrated material is used for the preparation of culture media. The manufacturers instructions relating to the preparation of these products should be rigorously followed.

5.2.2.2 Water

The water shall be freshly glass distilled water and not demineralized water.

Sterilize in the autoclave (see 5.3.1).

NOTE 1 If the water is sterilized during the sterilization of the reagents, this is not necessary.

NOTE 2 If distilled water of adequate quality is not available, water for injectable preparation (European Pharmacopoeia) can be used.

5.2.2.3 M17-broth

For maintenance of bacterial host strain, propagation of bacteriophages and for formulation of phage diluent (see 5.2.2.6).

Phytone peptone (from soya meal)	5,00 g
Polypeptone peptone (from casein & animal tissue)	5,00 g
Beef extract powder	5,00 g
Yeast extract	2,50 g
D(+)-lactose	5,00 g
Ascorbic acid	0,50 g
Sodium-β-glycerophosphate	19,00 g
Magnesium sulfate , 7 H ₂ O	0,25 g
Water (see 5.2.2.2)	1 000 ml

Sterilize in the autoclave (see 5.3.1). After sterilization the pH of the medium shall be equivalent to 7,0 ± 0,2 when measured at 20 °C. When M17-broth is the diluent for neutralizer formulation (see 5.2.2.12 and Annex B), double

1) DSM 4262, DSM 10567 and DSM 4366 are the collection numbers of bacteriophage and bacterial strains supplied by the DSMZ (Deutsche Sammlung von Mikroorganismen und ZellKulturen). This information is given for the convenience of users of this standard and does not constitute an endorsement by CEN of the culture collection named. Corresponding strains supplied by other culture collections may be used if they can be shown to lead to the same results.

concentrated M17-broth shall be used for preparation (i.e., all reagents shall be added in double concentration to 1 000 ml water).

5.2.2.4 M17-agar (bottom agar)

Bottom agar for quantitative counting of lysis zones (plaques) obtained from single infective bacteriophage particles in the bacterial lawn of the host bacteria.

Add 15 g of agar to 1 000 ml of M17-broth (see 5.2.2.3). Dissolve the agar by boiling with constant stirring.

Sterilize in the autoclave (see 5.3.1). After sterilization the pH of the medium shall be equivalent to $7,0 \pm 0,2$ when measured at 20 °C. When the agar is cooled down to $47 \text{ °C} \pm 1 \text{ °C}$, add 10 ml of a sterile 1 mol/l CaCl_2 -stock solution (see 5.2.2.8). Mix gently and pour 15 ml to 18 ml of agar into Petri dishes (see 5.3.2.10).

5.2.2.5 Overlay agar (top agar, soft agar)

For counting bacteriophages: Dissolve 6,5 g agar in 1 000 ml M17-broth (see 5.2.2.3) and heat until boiling with constant stirring. Dispense the molten agar in test tubes (2,5 to 3 ml each).

Sterilize in the autoclave (see 5.3.1).

NOTE For achieving clear phage-derived lysis zones (plaques) in the lawn of host bacterial cells only well-defined agar should be used which is specified by the supplier for phage enumeration by the overlay technique (see 5.5.2 and 5.5.3).

5.2.2.6 Phage diluent (on basis of ¼ strength Ringer's solution)

For preparing dilution series for titration of phage (counting of phage-derived lysis zones):

— 1/4-strength Ringer's solution:

— sodium chloride	2,250 g;
— potassium chloride	0,105 g;
— calcium chloride, anhydrous	0,06 g;
— sodium hydrogen carbonate	0,050 g;
— water (see 5.2.2.2)	to 1 000 ml.

Add 10 ml M17-broth (see 5.2.2.3) to 90 ml of 1/4-strength Ringer's solution.

Sterilize in the autoclave (see 5.3.1). Before use, add 1 ml from an 1 mol/l CaCl_2 -stock solution (see 5.2.2.8) to 100 ml of the dilution broth.

NOTE Ringer's solution can be prepared from ready-to-use tablets according to the supplier's recommendations.

5.2.2.7 SM-buffer

For resuspension and storage of intact phage particles:

— Tris-HCl	2,4 g;
— NaCl	5,8 g;