

SVENSK STANDARD

SS-EN 15742:2009

Fastställt/Approved: 2009-03-09

Publicerad/Published: 2009-04-07

Utgåva/Edition: 1

Språk/Language: engelska/English

ICS: 65.120

Djurfoder – Bestämning av organoklor-pesticider och PCB:er med GC/ECD

Animal feeding stuffs – Determination of OC-pesticides and PCB's by GC/ECD



SWEDISH
STANDARDS
INSTITUTE

Hitta rätt produkt och ett leveranssätt som passar dig

Standarder

Genom att följa gällande standard både effektiviserar och säkrar du ditt arbete. Många standarder ingår dessutom ofta i paket.

Tjänster

Abonnemang är tjänsten där vi uppdaterar dig med aktuella standarder när förändringar sker på dem du valt att abonnera på.

På så sätt är du säker på att du alltid arbetar efter rätt utgåva.

e-nav är vår online-tjänst som ger dig och dina kollegor tillgång till standarder ni valt att abonnera på dygnet runt. Med e-nav kan samma standard användas av flera personer samtidigt.

Leveranssätt

Du väljer hur du vill ha dina standarder levererade. Vi kan erbjuda dig dem på papper och som pdf.

Andra produkter

Vi har böcker som underlättar arbetet att följa en standard. Med våra böcker får du ökad förståelse för hur standarder ska följas och vilka fördelar den ger dig i ditt arbete. Vi tar fram många egna publikationer och fungerar även som återförsäljare. Det gör att du hos oss kan hitta över 500 unika titlar. Vi har även tekniska rapporter, specifikationer och "workshop agreement".

Matriser är en översikt på standarder och handböcker som bör läsas tillsammans. De finns på sis.se och ger dig en bra bild över hur olika produkter hör ihop.

Standardiseringsprojekt

Du kan påverka innehållet i framtida standarder genom att delta i någon av SIS ca 400 Tekniska Kommittéer.

Find the right product and the type of delivery that suits you

Standards

By complying with current standards, you can make your work more efficient and ensure reliability. Also, several of the standards are often supplied in packages.

Services

Subscription is the service that keeps you up to date with current standards when changes occur in the ones you have chosen to subscribe to. This ensures that you are always working with the right edition.

e-nav is our online service that gives you and your colleagues access to the standards you subscribe to 24 hours a day. With e-nav, the same standards can be used by several people at once.

Type of delivery

You choose how you want your standards delivered. We can supply them both on paper and as PDF files.

Other products

We have books that facilitate standards compliance. They make it easier to understand how compliance works and how this benefits you in your operation. We produce many publications of our own, and also act as retailers. This means that we have more than 500 unique titles for you to choose from. We also have technical reports, specifications and workshop agreements.

Matrices, listed at sis.se, provide an overview of which publications belong together.

Standardisation project

You can influence the content of future standards by taking part in one or other of SIS's 400 or so Technical Committees.

Europastandarden EN 15742:2009 gäller som svensk standard. Detta dokument innehåller den officiella engelska versionen av EN 15742:2009.

The European Standard EN 15742:2009 has the status of a Swedish Standard. This document contains the official English version of EN 15742:2009.

© Copyright/Upphovsrätten till denna produkt tillhör SIS, Swedish Standards Institute, Stockholm, Sverige. Användningen av denna produkt regleras av slutanvändarlicensen som återfinns i denna produkt, se standardens sista sidor.

© Copyright SIS, Swedish Standards Institute, Stockholm, Sweden. All rights reserved. The use of this product is governed by the end-user licence for this product. You will find the licence in the end of this document.

Upplysningar om sakinnehållet i standarden lämnas av SIS, Swedish Standards Institute, telefon 08-555 520 00. Standarder kan beställas hos SIS Förlag AB som även lämnar allmänna upplysningar om svensk och utländsk standard.

Information about the content of the standard is available from the Swedish Standards Institute (SIS), tel +46 8 555 520 00. Standards may be ordered from SIS Förlag AB, who can also provide general information about Swedish and foreign standards.

SIS Förlag AB, SE 118 80 Stockholm, Sweden. Tel: +46 8 555 523 10. Fax: +46 8 555 523 11.
E-mail: sis.sales@sis.se Internet: www.sis.se

EUROPEAN STANDARD

EN 15742

NORME EUROPÉENNE

EUROPÄISCHE NORM

February 2009

ICS 65.120

English Version

Animal feeding stuffs - Determination of OC-pesticides and PCB's by GC/ECD

Aliments des animaux - Détermination des pesticides organochlorés (OC) et des polychlorobiphényles (PCB) par GC/ECD

Futtermittel - Bestimmung der OC-Pestizide und PCB's mittels GC/ECD-Verfahren

This European Standard was approved by CEN on 24 January 2009.

CEN members are bound to comply with the CEN/CENELEC Internal Regulations which stipulate the conditions for giving this European Standard the status of a national standard without any alteration. Up-to-date lists and bibliographical references concerning such national standards may be obtained on application to the CEN Management Centre or to any CEN member.

This European Standard exists in three official versions (English, French, German). A version in any other language made by translation under the responsibility of a CEN member into its own language and notified to the CEN Management Centre has the same status as the official versions.

CEN members are the national standards bodies of Austria, Belgium, Bulgaria, Cyprus, Czech Republic, Denmark, Estonia, Finland, France, Germany, Greece, Hungary, Iceland, Ireland, Italy, Latvia, Lithuania, Luxembourg, Malta, Netherlands, Norway, Poland, Portugal, Romania, Slovakia, Slovenia, Spain, Sweden, Switzerland and United Kingdom.



EUROPEAN COMMITTEE FOR STANDARDIZATION
COMITÉ EUROPÉEN DE NORMALISATION
EUROPÄISCHES KOMITEE FÜR NORMUNG

Management Centre: Avenue Marnix 17, B-1000 Brussels

Contents

Page

Foreword	3
1 Scope	4
2 Normative references	4
3 Terms and definitions	4
4 Principle	5
5 Reagents and materials	5
6 Apparatus	12
7 Sampling	13
8 Preparation of test sample	13
9 Procedure	13
10 Calculation and expression of results	15
11 Precision	18
12 Test report	21
13 Important considerations	21
Annex A (informative) Results of interlaboratory tests	22
Bibliography	23

Foreword

This document (EN 15742:2009) has been prepared by Technical Committee CEN/TC 327 “Animal feeding stuffs”, the secretariat of which is held by NEN.

This European Standard shall be given the status of a national standard, either by publication of an identical text or by endorsement, at the latest by August 2009, and conflicting national standards shall be withdrawn at the latest by August 2009.

Attention is drawn to the possibility that some of the elements of this document may be the subject of patent rights. CEN [and/or CENELEC] shall not be held responsible for identifying any or all such patent rights.

This document has been prepared under a mandate given to CEN by the European Commission and the European Free Trade Association, and supports essential requirements of EC Directive(s).

According to the CEN/CENELEC Internal Regulations, the national standards organizations of the following countries are bound to implement this European Standard: Austria, Belgium, Bulgaria, Cyprus, Czech Republic, Denmark, Estonia, Finland, France, Germany, Greece, Hungary, Iceland, Ireland, Italy, Latvia, Lithuania, Luxembourg, Malta, Netherlands, Norway, Poland, Portugal, Romania, Slovakia, Slovenia, Spain, Sweden, Switzerland and the United Kingdom.

1 Scope

This European Standard specifies a gas chromatographic method with electron capture detection (ECD) for the determination of organochlorine pesticides (OC's) and polychlorinated biphenyls (PCBs) in animal feeding stuffs.

The method is applicable to animal feeding stuffs with a water content up to about 20 wt% and oil/fatty samples containing residues of one or more of the following OC's, PCBs, toxaphene and some of their isomers and degradation products:

- Aldrin;
- Dieldrin;
- Chlorocamphene (Toxaphene);
- Chlordane (= sum of Chlordane isomers and Oxychlordane);
- DDT (= sum of isomers op'-DDT, pp'-DDT, pp'-TDE (pp'-DDD), and pp'-DDE);
- Endosulfan (sum of α -/ β -isomers and Endosulfan-sulphate);
- Endrin;
- Heptachlor (= sum of Heptachlor and β -Heptachlorepoxyde);
- Hexachlorobenzene (HCB);
- Hexachlorocyclohexane isomers α -HCH (α -BHC), β -HCH (β -BHC), γ -HCH (γ -BHC or lindane);
- PCB 28, 52, 101, 138, 153 and 180 ("Indicator PCBs") and PCB 198, 209.

The limit of quantification for the mentioned organochlorine pesticides and PCBs is 5 ng/g in general. However, 10 ng/g applies for Heptachlor, Aldrin, Endrin, Dieldrin, and Endosulfan (α -, β - and sulphate). Individual laboratories are responsible to ensure that the equipment they used will achieve these limits of quantifications.

2 Normative references

The following referenced documents are indispensable for the application of this document. For dated references, only the edition cited applies. For undated references, the latest edition of the referenced document (including any amendments) applies.

ISO 6498, *Animal feeding stuffs – Preparation of test samples*

3 Terms and definitions

For the purpose of this document, the following terms and definitions apply.

3.1

limit of detection

smallest measured content, from which it is possible to deduce the presence of the analyte with reasonable statistical certainty

NOTE The limit of detection is numerically equal to three times the standard deviation of the mean of blank determinations ($n > 10$).

3.2

limit of quantification

lowest content of the analyte which can be measured with reasonable statistical certainty

NOTE If both accuracy and precision are constant over a concentration range around the limit of detection, then the limit of quantification is numerically equal to 6 times the standard deviation of the mean of blank determinations ($n > 10$).

3.3

feed additives

substances are feed additives when they comply with the definition of feed additives given in the Regulation 1831/2003

4 Principle

A test portion of animal feeding stuff is fortified with internal standard (PCB 198), and is extracted with ethylacetate. The extract is concentrated and subsequently purified by:

- Gel permeation chromatography (GPC), with cyclohexane/ethylacetate as eluting solvent
- chromatography on partially deactivated silica gel.

The collected fraction containing the compounds of interest is concentrated and re-dissolved in a solution containing another internal standard (PCB 209) as a reference standard. After cleanup the analytes are measured using GC-ECD. Identification is done on the basis of comparing retention times on capillary columns of different polarity. Quantification is done using the internal standard method.

5 Reagents and materials

Use only reagents of recognized analytical grade and with a purity suitable for OC and PCB residue analysis. Check the purity of the reagents by performing a blank test under the same conditions as used in the method. The chromatogram should not show any interfering impurity at the retention time of compounds of interest.

WARNING — The use of this European Standard can involve hazardous materials, operations and equipment. This standard does not purport to address all the safety problems associated with its use. It is the responsibility of the user of this European Standard to establish appropriate safety and health practices and determine the applicability of regulatory limitations prior to use.

5.1 Cyclohexane

5.2 Ethylacetate

5.3 Hexane

5.4 Dichloromethane

5.5 Iso-octane

5.6 Toluene

5.7 Hexane/toluene = 3+7, parts by volume

Mix 30 ml of hexan (5.3) with 70 ml of toluene (5.6) thoroughly. Store at room temperature in a tightly closed glass bottle.

5.8 Sodium Sulphate, anhydrous

Heated at 160°C during at least 24 h.

5.9 Ethylacetate/Cyclohexane = 1+1, parts by volume

Mix 500 ml of ethylacetate (5.2) with 500 ml of cyclohexane (5.1) thoroughly. Store at room temperature in a tightly closed glass bottle.

5.10 Silica gel, deactivated with 3,5% water

Heat silica gel 60 (63µm to 200µm = 70 mesh to 230 mesh), at 130°C for at least 5 h, allow to cool in a desiccator, and store in a tightly stoppered container in the desiccator. To 96,5 g dried silica gel in a 300 ml Erlenmeyer flask with a ground joint, add 3,5 ml water dropwise from a burette, with continuous swirling. Immediately stopper the flask with a ground stopper and shake vigorously for 5 min until all lumps have disappeared. Next shake for 2 h on a mechanical shaker, and then store in a tightly stoppered container. Deactivated silica gel is tenable during approximately 2 weeks if carefully stored.

5.11 Internal standard (PCB 198)

5.12 Internal Standard (PCB 209)

5.13 OC-pesticide reference standards

Each with a purity not less than 99%.

Aldrin

((1R,4S,4aS,5S,8R,8aR)-1,2,3,4,10,10-hexachloro-1,4,4a,5,8,8a-hexahydro-1,4:5,8-dimethanonaphthalene)
CAS Number: 309-00-2

Dieldrin

((1R,4S,4aS,5R,6R,7S,8S,8aR)-1,2,3,4,10,10-hexachloro-1,4,4a,5,6,7,8,8a-octahydro-6,7epoxy-1,4:5,8-dimethanonaphthalene)
CAS Number: 60-57-1

Chlordane

(1,2,4,5,6,7,8,8-octachloro-2,3,3a,4,7,7a-hexahydro-4,7-ethano-1*H*-indene); α and β isomer
CAS Numbers: 5103-71-9 and 5103-74-2

Oxychlordane

(4,7-Methanoindan, 1,2,4,5,6,7,8,8-octachloro-2,3-epoxy-3a,4,7,7a-tetrahydro-, exo,endo-)
CAS Number: 27304-13-8

op'-DDT

[o,p'-(1,1,1-trichloro-2,2-bis(4-chlorophenyl)ethane)]
CAS Number: 789-02-6

pp'-DDT

[p,p'-(1,1,1-trichloro-2,2-bis(4-chlorophenyl) ethane)]

CAS Number: 50-29-3

pp'-TDE

(p,p'-DDD) [p,p'-1,1-dichloro-2,2-bis(4-chlorophenyl) ethane]

CAS Number: 72-54-8

pp'-DDE

[p,p'-(1,1-dichloro-2,2-bis(4-chlorophenyl) ethylene)]

CAS Number: 72-55-9

Endosulfan

(6,7,8,9,10,10-hexachloro-1,5,5a,6,9,9a-hexahydro-6,9-methano-2,4,3-benzodioxathiepin 3-oxide)

two stereoisomers, α , (I), CAS Number: 959-98-8 and β , (II), CAS Number: 33213-65-9.

Endosulfan-sulphate;

CAS Number: 1031-07-8

Endrin

[(1R,4S,4aS,5S,6S,7R,8R,8aR)-1,2,3,4,10,10-hexachloro-1,4,4a,5,6,7,8,8a-octahydro-6,7-epoxy-1,4:5,8-dimethanonaphthalene]

CAS Number: 72-20-8

Heptachlor

(1,4,5,6,7,8,8-heptachloro-3a,4,7,7a-tetrahydro-4,7-methanoindene)

CAS Number: 76-44-8

β -Heptachlorepoxyde

(1,4,5,6,7,8,8-heptachloro-3a,4,7,7a-tetrahydro-4,7-methanoindene(exo))

CAS Number: 1024-57-3

HCB

(hexachlorobenzene)

CAS Number: 118-74-1

α -HCH (α -BHC)

(α -1,2,3,4,5,6-hexachlorocyclohexane)

CAS Number: 319-84-6

β -HCH (β -BHC)

(β -1,2,3,4,5,6-hexachlorocyclohexane)

CAS Number: 319-85-7

γ -HCH (γ -BHC; lindane)

(1,2,3,4,5,6-hexachlorocyclohexane)

CAS Number: 58-89-9

Or a Certified Mixture at a concentration of 10 $\mu\text{g/ml}$.

5.14 PCBs reference standards

Each with a purity not less than 99%.

PCB 28 (2,4,4' trichlorobiphenyl) ; CAS Number: 7012-37-5

PCB 52 (2,2',5,5' tetrachlorobiphenyl) ; CAS Number: 35693-99-3