

# SVENSK STANDARD

## SS-EN 15652:2009

Fastställt/Approved: 2009-06-05  
Publicerad/Published: 2009-08-18  
Utgåva/Edition: 1  
Språk/Language: engelska/English  
ICS: 67.050

---

### Foodstuffs – Determination of niacin by HPLC



# Hitta rätt produkt och ett leveranssätt som passar dig

## Standarder

Genom att följa gällande standard både effektiviserar och säkrar du ditt arbete. Många standarder ingår dessutom ofta i paket.

## Tjänster

Abonnemang är tjänsten där vi uppdaterar dig med aktuella standarder när förändringar sker på dem du valt att abonnera på.

På så sätt är du säker på att du alltid arbetar efter rätt utgåva.

e-nav är vår online-tjänst som ger dig och dina kollegor tillgång till standarder ni valt att abonnera på dygnet runt. Med e-nav kan samma standard användas av flera personer samtidigt.

## Leveranssätt

Du väljer hur du vill ha dina standarder levererade. Vi kan erbjuda dig dem på papper och som pdf.

## Andra produkter

Vi har böcker som underlättar arbetet att följa en standard. Med våra böcker får du ökad förståelse för hur standarder ska följas och vilka fördelar den ger dig i ditt arbete. Vi tar fram många egna publikationer och fungerar även som återförsäljare. Det gör att du hos oss kan hitta över 500 unika titlar. Vi har även tekniska rapporter, specifikationer och "workshop agreement".

Matriser är en översikt på standarder och handböcker som bör läsas tillsammans. De finns på sis.se och ger dig en bra bild över hur olika produkter hör ihop.

## Standardiseringsprojekt

Du kan påverka innehållet i framtida standarder genom att delta i någon av SIS ca 400 Tekniska Kommittéer.

# Find the right product and the type of delivery that suits you

## Standards

By complying with current standards, you can make your work more efficient and ensure reliability. Also, several of the standards are often supplied in packages.

## Services

Subscription is the service that keeps you up to date with current standards when changes occur in the ones you have chosen to subscribe to. This ensures that you are always working with the right edition.

e-nav is our online service that gives you and your colleagues access to the standards you subscribe to 24 hours a day. With e-nav, the same standards can be used by several people at once.

## Type of delivery

You choose how you want your standards delivered. We can supply them both on paper and as PDF files.

## Other products

We have books that facilitate standards compliance. They make it easier to understand how compliance works and how this benefits you in your operation. We produce many publications of our own, and also act as retailers. This means that we have more than 500 unique titles for you to choose from. We also have technical reports, specifications and workshop agreements.

Matrices, listed at sis.se, provide an overview of which publications belong together.

## Standardisation project

You can influence the content of future standards by taking part in one or other of SIS's 400 or so Technical Committees.

Europastandarden EN 15652:2009 gäller som svensk standard. Detta dokument innehåller den officiella engelska versionen av EN 15652:2009.

The European Standard EN 15652:2009 has the status of a Swedish Standard. This document contains the official English version of EN 15652:2009.

© Copyright/Upphovsrätten till denna produkt tillhör SIS, Swedish Standards Institute, Stockholm, Sverige. Användningen av denna produkt regleras av slutanvändarlicensen som återfinns i denna produkt, se standardens sista sidor.

© Copyright SIS, Swedish Standards Institute, Stockholm, Sweden. All rights reserved. The use of this product is governed by the end-user licence for this product. You will find the licence in the end of this document.

Upplysningar om sakinnehållet i standarden lämnas av SIS, Swedish Standards Institute, telefon 08-555 520 00.

Standarder kan beställas hos SIS Förlag AB som även lämnar allmänna upplysningar om svensk och utländsk standard.

Information about the content of the standard is available from the Swedish Standards Institute (SIS), tel +46 8 555 520 00.

Standards may be ordered from SIS Förlag AB, who can also provide general information about Swedish and foreign standards.

SIS Förlag AB, SE 118 80 Stockholm, Sweden. Tel: +46 8 555 523 10. Fax: +46 8 555 523 11.

E-mail: [sis.sales@sis.se](mailto:sis.sales@sis.se) Internet: [www.sis.se](http://www.sis.se)



EUROPEAN STANDARD  
NORME EUROPÉENNE  
EUROPÄISCHE NORM

**EN 15652**

May 2009

---

ICS 67.050

English Version

## Foodstuffs - Determination of niacin by HPLC

Produits alimentaires - Dosage de la niacine par CLHP

Lebensmittel - Bestimmung von Niacin mit HPLC

This European Standard was approved by CEN on 23 April 2009.

CEN members are bound to comply with the CEN/CENELEC Internal Regulations which stipulate the conditions for giving this European Standard the status of a national standard without any alteration. Up-to-date lists and bibliographical references concerning such national standards may be obtained on application to the CEN Management Centre or to any CEN member.

This European Standard exists in three official versions (English, French, German). A version in any other language made by translation under the responsibility of a CEN member into its own language and notified to the CEN Management Centre has the same status as the official versions.

CEN members are the national standards bodies of Austria, Belgium, Bulgaria, Cyprus, Czech Republic, Denmark, Estonia, Finland, France, Germany, Greece, Hungary, Iceland, Ireland, Italy, Latvia, Lithuania, Luxembourg, Malta, Netherlands, Norway, Poland, Portugal, Romania, Slovakia, Slovenia, Spain, Sweden, Switzerland and United Kingdom.



EUROPEAN COMMITTEE FOR STANDARDIZATION  
COMITÉ EUROPÉEN DE NORMALISATION  
EUROPÄISCHES KOMITEE FÜR NORMUNG

**Management Centre: Avenue Marnix 17, B-1000 Brussels**

<b>Contents</b>		Page
<b>Foreword</b> .....		<b>3</b>
<b>1</b>	<b>Scope</b> .....	<b>4</b>
<b>2</b>	<b>Normative references</b> .....	<b>4</b>
<b>3</b>	<b>Principle</b> .....	<b>4</b>
<b>4</b>	<b>Reagents</b> .....	<b>4</b>
<b>5</b>	<b>Apparatus</b> .....	<b>7</b>
<b>6</b>	<b>Procedure</b> .....	<b>10</b>
<b>7</b>	<b>Calculation</b> .....	<b>11</b>
<b>8</b>	<b>Precision</b> .....	<b>12</b>
<b>9</b>	<b>Test report</b> .....	<b>14</b>
<b>Annex A</b> (informative) <b>Typical chromatogram</b> .....		<b>15</b>
<b>Annex B</b> (informative) <b>Precision data for acid-, enzymatic- and acid/alkaline hydrolysis</b> .....		<b>16</b>
<b>Annex C</b> (informative) <b>Comparison between three different ways of hydrolysis</b> .....		<b>19</b>
<b>Bibliography</b> .....		<b>21</b>

## Foreword

This document (EN 15652:2009) has been prepared by Technical Committee CEN/TC 275 “Food analysis - Horizontal methods”, the secretariat of which is held by DIN.

This European Standard shall be given the status of a national standard, either by publication of an identical text or by endorsement, at the latest by November 2009, and conflicting national standards shall be withdrawn at the latest by November 2009.

Attention is drawn to the possibility that some of the elements of this document may be the subject of patent rights. CEN [and/or CENELEC] shall not be held responsible for identifying any or all such patent rights.

**WARNING —The use of this standard may involve hazardous materials, operations and equipment. This standard does not purport to address all the safety problems associated with its use. It is the responsibility of the user of this standard to establish appropriate safety and health practices and determine the applicability of regulatory limitations prior to use.**

According to the CEN/CENELEC Internal Regulations, the national standards organizations of the following countries are bound to implement this European Standard: Austria, Belgium, Bulgaria, Cyprus, Czech Republic, Denmark, Estonia, Finland, France, Germany, Greece, Hungary, Iceland, Ireland, Italy, Latvia, Lithuania, Luxembourg, Malta, Netherlands, Norway, Poland, Portugal, Romania, Slovakia, Slovenia, Spain, Sweden, Switzerland and the United Kingdom.

## 1 Scope

This European Standard specifies a method for the determination of the mass fraction of niacin in foodstuffs by high performance liquid chromatography (HPLC) by three different ways of hydrolysis, acid hydrolysis (A), enzymatic hydrolysis (B) or acid/alkaline hydrolysis (C).

The method has been validated in interlaboratory tests on fortified and non-fortified samples such as breakfast cereal powder, chocolate cereals, cooked ham, green peas, lyophilized green peas with ham, lyophilized soup, nutritive orange juice, milk powder and wheat flour, at levels from 0,5 mg/100 g to 24 mg/100 g. For further information on the validation data, see Annex B.

A and B give similar results for niacin. In options A and B niacin is calculated as the sum of nicotinamide and nicotinic acid, and expressed as nicotinic acid [1]. Option C gives higher results than A and B for niacin with non-supplemented cereals, but similar results for other products. In option C, niacin is calculated and expressed as nicotinic acid after transformation of nicotinamide into nicotinic acid [2].

Option A is faster and cheaper than B and C.

Option B is used if an exact quantification of nicotinamide and nicotinic acid is needed. This cannot be done with option A, because there is a slight transformation of nicotinamide into nicotinic acid during the acid hydrolysis.

Option C quantifies total niacin. The alkaline hydrolysis is able to liberate other forms giving higher results for niacin, which in some foods such as maize and cereals are not normally biologically available, see [3], [4] and [5].

Information on a comparison between the three different ways of hydrolysis is given in Annex C.

## 2 Normative references

The following referenced documents are indispensable for the application of this document. For dated references, only the edition cited applies. For undated references, the latest edition of the referenced document (including any amendments) applies.

EN ISO 3696:1995, *Water for analytical laboratory use – Specification and test methods (ISO 3696:1987)*

## 3 Principle

Niacin vitamers are extracted from food by an acid (option A), an enzymatic (option B) or an acid/alkaline (option C) treatment and quantified by HPLC with a fluorimetric detection after a post-column derivatization with UV irradiation, see [1] and [2]. For option A and option B, niacin is determined as the sum of nicotinamide and nicotinic acid. Niacin is expressed as nicotinic acid after correction of the molecular weights. For option C, niacin is determined and expressed as nicotinic acid. The alkaline treatment transforms all nicotinamide into nicotinic acid.

## 4 Reagents

### 4.1 General

During the analysis, unless otherwise stated, use only reagents of recognized analytical grade and water of at least grade 1 according to EN ISO 3696:1995.



## 4.2 Chemicals and solutions

- 4.2.1 Sodium acetate**, mass fraction  $w(\text{C}_2\text{H}_3\text{NaO}_2) \geq 99 \%$
- 4.2.2 Potassium monohydrogen phosphate**,  $w(\text{K}_2\text{HPO}_4) \geq 99,5 \%$
- 4.2.3 Potassium dihydrogen phosphate**,  $w(\text{KH}_2\text{PO}_4) \geq 99,5 \%$
- 4.2.4 Non stabilized hydrogen peroxyde solution**,  $w(\text{H}_2\text{O}_2) = 30 \%$
- 4.2.5 Copper sulfate**,  $w(\text{Cu(II)SO}_4 \cdot 5\text{H}_2\text{O}) \geq 99 \%$
- 4.2.6 Acetic acid**,  $w(\text{CH}_3\text{COOH}) \geq 99,8 \%$
- 4.2.7 Concentrated hydrochloric acid solution (option A and C)**,  $w(\text{HCl}) = 37,0 \%$
- 4.2.8 NADase from *Neurospora crassa* (option B)**, enzyme activity **0,55 U/mg of protein**.

Store below 0 °C.

NOTE For the interlaboratory study, NADase from *Neurospora crassa* from Sigma Chemicals with reference N 9629, lyophilised powder, 0,5 U/mg to 3,0 U/mg protein (biuret) has been used<sup>1</sup>.

- 4.2.9 Acetic acid solution**, substance concentration  $c(\text{CH}_3\text{COOH}) = 5 \text{ mol/l}$
- 4.2.10 Sodium acetate solution**,  $c(\text{C}_2\text{H}_3\text{NaO}_2) = 2,5 \text{ mol/l}$
- 4.2.11 Copper sulfate solution**,  $c(\text{Cu(II)SO}_4 \cdot 5\text{H}_2\text{O}) = 0,005 \text{ mol/l}$
- 4.2.12 Sodium acetate solution (option B)**,  $c(\text{C}_2\text{H}_3\text{NaO}_2) = 0,05 \text{ mol/l}$ , pH = 4,5

Dissolve 4,10 g of sodium acetate (4.2.1) in 900 ml of water. Adjust the solution to pH = 4,5 with acid acetic (4.2.6), and then dilute to 1000 ml with water.

- 4.2.13 Phosphate buffer solution (option B)**,  $c(\text{K}_2\text{HPO}_4) = 0,05 \text{ mol/l}$  and  $c(\text{KH}_2\text{PO}_4) = 0,05 \text{ mol/l}$ , pH = 6,8

Mix 1 part per volume of  $\text{K}_2\text{HPO}_4$  solution (0,05 mol/l) and 1 part per volume of  $\text{KH}_2\text{PO}_4$  solution (0,05 mol/l). Adjust pH to 6,8 with sodium acetate solution (4.2.10) if necessary.

- 4.2.14 NADase solution (option B)**

Dissolve 2,9 mg of NADase (4.2.8) in 5 ml of phosphate buffer (4.2.13). This solution is stable for 1 week at -18 °C.

- 4.2.15 Hydrochloric acid solution (options A and C)**,  $c(\text{HCl}) = 0,1 \text{ mol/l}$

---

<sup>1</sup> This information is given for the convenience of users of this European Standard and does not constitute an endorsement by CEN of the product named. Equivalent products may be used if they can be shown to lead to the same results.

#### 4.2.16 HPLC mobile phase

Dissolve 4,77 g of potassium dihydrogen phosphate (4.2.3) in 400 ml of water. Add 3,8 ml of hydrogen peroxide solution (4.2.4) and 0,5 ml of copper sulfate solution (4.2.11). Dilute to 500 ml. The pH is about 4,5. Filter through a membrane filter (5.7). This solution is stable for one day.

#### 4.2.17 Sodium hydroxide (option C), $w(\text{NaOH}) \geq 99\%$

#### 4.2.18 Sodium hydroxide solution (option C), $c(\text{NaOH}) = 5\text{ mol/l}$

Dissolve 20 g of sodium hydroxide (4.2.17) in 80 ml of water. After cooling dilute to 100 ml.

#### 4.2.19 Hydrochloric acid solution (option C), $w(\text{HCl}) = 3,7\%$

Dilute 5 ml of the concentrated hydrochloric acid solution (4.2.7) to 50 ml with water.

### 4.3 Standard substances

#### 4.3.1 General

Nicotinic acid and nicotinamide can be obtained from various suppliers. The purity may vary and it is therefore necessary to determine the concentration of the calibration solution by a spectrometric determination (see 4.4.3).

#### 4.3.2 Nicotinic acid, $w(\text{C}_6\text{H}_5\text{NO}_2) \geq 99,5\%$

#### 4.3.3 Nicotinamide (options A and B), $w(\text{C}_6\text{H}_6\text{N}_2\text{O}) \geq 99,5\%$

### 4.4 Stock solutions

#### 4.4.1 Nicotinic acid stock solution, mass concentration $\rho(\text{C}_6\text{H}_5\text{NO}_2) = 1\text{ mg/ml}$

Dissolve an amount of the nicotinic acid standard substance (4.3.2), e.g. approximately 100 mg to the nearest 1 mg in 100 ml of water. This solution is stable for 1 week at  $-18\text{ }^\circ\text{C}$ .

#### 4.4.2 Nicotinamide stock solution, (options A and B), $\rho(\text{C}_6\text{H}_6\text{N}_2\text{O}) = 1\text{ mg/ml}$

Dissolve an amount of the nicotinamide standard substance (4.3.3), e.g. approximately 100 mg to the nearest 1 mg in 100 ml of water. This solution is stable for 1 week at  $-18\text{ }^\circ\text{C}$ .

### 4.4.3 Concentration tests

#### 4.4.3.1 Nicotinic acid solution, $\rho(\text{C}_6\text{H}_5\text{NO}_2) = 1\text{ mg/ml}$

Dilute 1 ml of the nicotinic acid stock solution (4.4.1) in 100 ml of hydrochloric acid solution (4.2.15) and measure the absorbance at 260 nm in a 1 cm cell using a UV spectrometer (5.2) against hydrochloric acid solution (4.2.15) as reference. Calculate the mass concentration,  $\rho$ , in milligram per millilitre of the stock solution, using Equation (1):

$$\rho = \frac{A_{260} \times 1000}{420} \quad (1)$$

where

$A_{260}$  is the absorbance value of the solution at 260 nm;

420 is the  $E_{1\text{cm}}^{1\%}$  value for nicotinic acid in 0,1 mol/l HCl, see [6].

#### 4.4.3.2 Nicotinamide solution, $\rho(\text{C}_6\text{H}_6\text{N}_2\text{O}) = 1 \text{ mg/ml}$

Dilute 1 ml of the nicotinamide stock solution (4.4.2) in 100 ml of hydrochloric acid solution (4.2.15) and measure the absorbance at 260 nm in a 1 cm cell using a UV spectrometer (5.2) against hydrochloric acid solution (4.2.15) as reference. Calculate the mass concentration,  $\rho$ , in milligram per millilitre of the stock solution, using Equation (2):

$$\rho = \frac{A_{260} \times 1000}{410} \quad (2)$$

where

$A_{260}$  is the absorbance value of the solution at 260 nm;

410 is the  $E_{1\text{cm}}^{1\%}$  value for nicotinamide in 0,1 mol/l HCl, see [6].

#### 4.5 Nicotinic acid and nicotinamide standard solutions, $\rho(\text{C}_6\text{H}_5\text{NO}_2) = \rho(\text{C}_6\text{H}_6\text{N}_2\text{O}) = 0,05 \text{ }\mu\text{g/ml}$ to 5 $\mu\text{g/ml}$

Prepare e.g. a first solution with 1 ml of each stock solution (4.4.1) or (4.4.2) in 100 ml of water. From this solution prepare four standard solutions (0,5 ml, 2,5 ml, 10 ml and 50 ml) in 100 ml of water. These solutions are stable for one day at room temperature.

## 5 Apparatus

### 5.1 General

Usual laboratory apparatus and glassware, and the following.

### 5.2 UV spectrometer

Capable of measurement of absorbance at defined wavelengths

### 5.3 Oven, capable of maintaining a temperature of 37 °C

### 5.4 Autoclave, capable of maintaining a temperature of 120 °C

### 5.5 HPLC system

Consisting of a pump, sample injecting device, fluorescence detector with excitation and emission wavelengths set at 322 nm and 380 nm, and an evaluation system such as an integrator.