

SVENSK STANDARD

SS-EN ISO 17380:2013

Fastställt/Approved: 2013-06-23
Publicerad/Published: 2013-06-26
Utgåva/Edition: 1
Språk/Language: engelska/English
ICS: 13.080.10

Markundersökningar – Bestämning av total cyanid och lättillgänglig cyanid – Kontinuerligt flöde analysmetod (ISO 17380:2013)

Soil quality – Determination of total cyanide and easily liberatable cyanide – Continuous-flow analysis method (ISO 17380:2013)

This preview is downloaded from www.sis.se. Buy the entire standard via <https://www.sis.se/std-98545>

Standarder får världen att fungera

SIS (Swedish Standards Institute) är en fristående ideell förening med medlemmar från både privat och offentlig sektor. Vi är en del av det europeiska och globala nätverk som utarbetar internationella standarder. Standarder är dokumenterad kunskap utvecklad av framstående aktörer inom industri, näringsliv och samhälle och befrämjar handel över gränser, bidrar till att processer och produkter blir säkrare samt effektiviserar din verksamhet.

Delta och påverka

Som medlem i SIS har du möjlighet att påverka framtida standarder inom ditt område på nationell, europeisk och global nivå. Du får samtidigt tillgång till tidig information om utvecklingen inom din bransch.

Ta del av det färdiga arbetet

Vi erbjuder våra kunder allt som rör standarder och deras tillämpning. Hos oss kan du köpa alla publikationer du behöver – allt från enskilda standarder, tekniska rapporter och standardpaket till handböcker och onlinetjänster. Genom vår webbtjänst e-nav får du tillgång till ett lättnavigerat bibliotek där alla standarder som är aktuella för ditt företag finns tillgängliga. Standarder och handböcker är källor till kunskap. Vi säljer dem.

Utveckla din kompetens och lyckas bättre i ditt arbete

Hos SIS kan du gå öppna eller företagsinterna utbildningar kring innehåll och tillämpning av standarder. Genom vår närhet till den internationella utvecklingen och ISO får du rätt kunskap i rätt tid, direkt från källan. Med vår kunskap om standarders möjligheter hjälper vi våra kunder att skapa verklig nytta och lönsamhet i sina verksamheter.

Vill du veta mer om SIS eller hur standarder kan effektivisera din verksamhet är du välkommen in på www.sis.se eller ta kontakt med oss på tel 08-555 523 00.



Standards make the world go round

SIS (Swedish Standards Institute) is an independent non-profit organisation with members from both the private and public sectors. We are part of the European and global network that draws up international standards. Standards consist of documented knowledge developed by prominent actors within the industry, business world and society. They promote cross-border trade, they help to make processes and products safer and they streamline your organisation.

Take part and have influence

As a member of SIS you will have the possibility to participate in standardization activities on national, European and global level. The membership in SIS will give you the opportunity to influence future standards and gain access to early stage information about developments within your field.

Get to know the finished work

We offer our customers everything in connection with standards and their application. You can purchase all the publications you need from us - everything from individual standards, technical reports and standard packages through to manuals and online services. Our web service e-nav gives you access to an easy-to-navigate library where all standards that are relevant to your company are available. Standards and manuals are sources of knowledge. We sell them.

Increase understanding and improve perception

With SIS you can undergo either shared or in-house training in the content and application of standards. Thanks to our proximity to international development and ISO you receive the right knowledge at the right time, direct from the source. With our knowledge about the potential of standards, we assist our customers in creating tangible benefit and profitability in their organisations.

If you want to know more about SIS, or how standards can streamline your organisation, please visit www.sis.se or contact us on phone +46 (0)8-555 523 00



Europastandarden EN ISO 17380:2013 gäller som svensk standard. Detta dokument innehåller den officiella engelska versionen av EN ISO 17380:2013.

The European Standard EN ISO 17380:2013 has the status of a Swedish Standard. This document contains the official version of EN ISO 17380:2013.

© Copyright/Upphovsrätten till denna produkt tillhör SIS, Swedish Standards Institute, Stockholm, Sverige. Användningen av denna produkt regleras av slutanvändarlicensen som återfinns i denna produkt, se standardens sista sidor.

© Copyright SIS, Swedish Standards Institute, Stockholm, Sweden. All rights reserved. The use of this product is governed by the end-user licence for this product. You will find the licence in the end of this document.

Upplysningar om sakinnehållet i standarden lämnas av SIS, Swedish Standards Institute, telefon 08-555 520 00. Standarder kan beställas hos SIS Förlag AB som även lämnar allmänna upplysningar om svensk och utländsk standard.

Information about the content of the standard is available from the Swedish Standards Institute (SIS), telephone +46 8 555 520 00. Standards may be ordered from SIS Förlag AB, who can also provide general information about Swedish and foreign standards.

Denna standard är framtagen av kommittén för Karaktärisering av avfall, mark och slam, SIS/TK 535.

Har du synpunkter på innehållet i den här standarden, vill du delta i ett kommande revideringsarbete eller vara med och ta fram andra standarder inom området? Gå in på www.sis.se - där hittar du mer information.

EUROPEAN STANDARD

EN ISO 17380

NORME EUROPÉENNE

EUROPÄISCHE NORM

June 2013

ICS 13.080.10

English Version

Soil quality - Determination of total cyanide and easily liberatable cyanide - Continuous-flow analysis method (ISO 17380:2013)

Qualité du sol - Détermination des cyanures totaux et des cyanures aisément libérables - Méthode d'analyse en flux continu (ISO 17380:2013)

Bodenbeschaffenheit - Bestimmung des Gehalts an gesamtem Cyanid und leicht freisetzbarem Cyanid - Verfahren mittels kontinuierlicher Durchflußanalyse (ISO 17380:2013)

This European Standard was approved by CEN on 8 June 2013.

CEN members are bound to comply with the CEN/CENELEC Internal Regulations which stipulate the conditions for giving this European Standard the status of a national standard without any alteration. Up-to-date lists and bibliographical references concerning such national standards may be obtained on application to the CEN-CENELEC Management Centre or to any CEN member.

This European Standard exists in three official versions (English, French, German). A version in any other language made by translation under the responsibility of a CEN member into its own language and notified to the CEN-CENELEC Management Centre has the same status as the official versions.

CEN members are the national standards bodies of Austria, Belgium, Bulgaria, Croatia, Cyprus, Czech Republic, Denmark, Estonia, Finland, Former Yugoslav Republic of Macedonia, France, Germany, Greece, Hungary, Iceland, Ireland, Italy, Latvia, Lithuania, Luxembourg, Malta, Netherlands, Norway, Poland, Portugal, Romania, Slovakia, Slovenia, Spain, Sweden, Switzerland, Turkey and United Kingdom.



EUROPEAN COMMITTEE FOR STANDARDIZATION
COMITÉ EUROPÉEN DE NORMALISATION
EUROPÄISCHES KOMITEE FÜR NORMUNG

Management Centre: Avenue Marnix 17, B-1000 Brussels

Contents		Page
Foreword		iv
Introduction		v
1 Scope		1
2 Normative references		1
3 Terms and definitions		1
4 Principle		2
4.1 Pretreatment of soil samples.....		2
4.2 Determination of total cyanide content.....		2
4.3 Determination of easily-liberatable cyanide content.....		2
5 Reagents		2
5.1 General reagents.....		2
5.2 Reagents for determination of cyanide.....		3
5.3 Standard solutions for cyanide determination.....		4
6 Apparatus		6
6.1 Standard laboratory glassware and laboratory equipment.....		6
6.2 Continuous flow analysis system.....		6
6.3 Shaking machine.....		7
7 Sampling, sample preparation and extraction		8
8 Working range		8
9 Procedure		8
9.1 General.....		8
9.2 Checking analysis system for correct function.....		8
9.3 Calibration graph.....		9
9.4 Measurement of samples.....		9
9.5 Cleaning of the analysis system.....		10
10 Calculation of the cyanide content in the sample		10
11 Precision		10
12 Test report		11
Annex A (informative) Recovery and repeatability		12
Annex B (informative) Information on continuous flow automatic distillation apparatus		14
Annex C (informative) Alternative distillation buffer (pH = 3,8)		15
Bibliography		16

Foreword

This document (EN ISO 17380:2013) has been prepared by Technical Committee ISO/TC 190 "Soil quality" in collaboration with Technical Committee CEN/TC 345 "Characterization of soils" the secretariat of which is held by NEN.

This European Standard shall be given the status of a national standard, either by publication of an identical text or by endorsement, at the latest by December 2013, and conflicting national standards shall be withdrawn at the latest by December 2013.

Attention is drawn to the possibility that some of the elements of this document may be the subject of patent rights. CEN [and/or CENELEC] shall not be held responsible for identifying any or all such patent rights.

According to the CEN-CENELEC Internal Regulations, the national standards organizations of the following countries are bound to implement this European Standard: Austria, Belgium, Bulgaria, Croatia, Cyprus, Czech Republic, Denmark, Estonia, Finland, Former Yugoslav Republic of Macedonia, France, Germany, Greece, Hungary, Iceland, Ireland, Italy, Latvia, Lithuania, Luxembourg, Malta, Netherlands, Norway, Poland, Portugal, Romania, Slovakia, Slovenia, Spain, Sweden, Switzerland, Turkey and the United Kingdom.

Endorsement notice

The text of ISO 17380:2013 has been approved by CEN as EN ISO 17380:2013 without any modification.

Introduction

Cyanides may be present in soil as cyanide ions and as complex cyanides. They can be determined as easily-liberatable cyanide or as total cyanide. Complex cyanide can be calculated by subtracting the easily-liberatable cyanide result from the total cyanide result. This International Standard specifies the determination of easily-liberatable cyanide, complex cyanides and total cyanide.

Methods using flow analysis automate wet chemical procedures and are particularly suitable for the processing of many analytes in water or soil extracts in large sample series at a high analysis frequency. The continuous flow analysis (CFA) method uses an automated dosage of the sample into a flow system (manifold) where the analytes in the sample react with the reagent solution on their way through the manifold. The sample preparation may be integrated in the manifold. The reaction product is measured in a photometric detector (e.g. flow cell photometer).^{[1],[2]}

In ISO 11262 a manual method for the photometric and volumetric determination of total cyanide in soil samples is described. It should be noted that the total cyanide results in soil samples as described in ISO 11262 may show slight differences from this International Standard. These differences are not considered to be very significant for this analysis. The easily-liberatable cyanide test has been removed from ISO 11262 because the validation data for this method were very poor.

For the analysis of cyanide in water ISO 14403-1^[6] and ISO 14403-2^[7] can be applied. The analytical procedure described in ISO 14403-2^[7] is identical to the one specified in this International Standard.

Soil quality — Determination of total cyanide and easily liberatable cyanide — Continuous-flow analysis method

WARNING — Cyanide solutions are highly toxic. Appropriate measures shall be taken to avoid ingestion. Care should be taken in the disposal of these solutions.

1 Scope

This International Standard specifies a method for the photometric determination of the total cyanide and easily-liberatable cyanide content in soil using automated distillation/continuous-flow analysis.

The International Standard applies to all types of soil with cyanide contents above 1 mg/kg on the basis of dry matter, expressed as cyanide ion.

NOTE Sulfide concentrations in the sample higher than 40 mg/kg dry matter cause interference. This effect can be recognized by the split peaks and as a slow decrease of the detector signal and can only be prevented by diluting the sample extract.

2 Normative references

The following documents, in whole or in part, are normatively referenced in this document and are indispensable for its application. For dated references, only the edition cited applies. For undated references, the latest edition of the referenced document (including any amendments) applies.

ISO 3696, *Water for analytical laboratory use — Specification and test methods*

ISO 9297, *Water quality — Determination of chloride — Silver nitrate titration with chromate indicator (Mohr's method)*

ISO 11262, *Soil quality — Determination of total cyanide*

ISO 11465, *Soil quality — Determination of dry matter and water content on a mass basis — Gravimetric method*

ISO 14507:2003, *Soil quality — Pretreatment of samples for determination of organic contaminants*

3 Terms and definitions

For the purposes of this document, the following terms and definitions apply.

3.1

total cyanide content

content of inorganic cyanide compounds consisting of the sum of the contents of easily-liberatable cyanide species and cyanide bound to metal cyanides, with the exception of thiocyanate ions and only partial breakdown of the cyanide bound in cobalt, gold, palladium and platinum cyanide complexes

3.2

easily-liberatable cyanide

content of inorganic cyanide compounds consisting of the sum of the contents of the easily-liberatable cyanide ions and the cyanide bound in simple metal cyanides (all expressed as CN) which are determined under the conditions of the method described in this International Standard

Note 1 to entry: The weakly complexed cyanide contained in tetracyanonickelate(II) and dicyanomercurate(II) is determined with the method for easily-liberatable cyanide. Up to 5 % of the strongly complexed cyanide in iron(III) hexacyanoferrate(II), hexacyanoferrate(III) and hexocyanoferrate(II) is determined with the method for easily-liberatable cyanide. Organic cyanide compounds (such as acetonitrile) are not determined.

Note 2 to entry: In connection with waste from the production of gold, "easily liberatable cyanide" is named "weak acid dissociable cyanide"^[8].

3.3 complex cyanide content

total cyanide content less than the easily-liberatable cyanide content

4 Principle

4.1 Pretreatment of soil samples

The soil sample is extracted with 2,5 mol/l sodium hydroxide solution for 16 h. The extract is diluted 100 times, after which analysis of total and/or easily-liberatable cyanide is performed.

4.2 Determination of total cyanide content

Complex bound cyanide, present in the diluted extract (0,025 mol/l NaOH), is decomposed in a continuously flowing stream after addition of a buffer solution, with an initial pH of 3,8, by the effect of UV light. An UV-B lamp and decomposition coil of borosilicate glass are used. UV light with a wavelength of less than 290 nm is absorbed by the glass, preventing the photolytic conversion of thiocyanate into cyanide. The hydrogen cyanide present at these conditions is separated by in-line distillation at a heating bath temperature of $125\text{ °C} \pm 2\text{ °C}$ and then determined photometrically. The photometric determination is based on the reaction of cyanide with chloramine-T resulting in the formation of cyanogen chloride. This reacts with pyridine-4-carboxylic acid and 1,3-dimethylbarbituric acid to give a purple colour. The absorbance at 606 nm is then measured to determine the cyanide content. When an interference filter is used, a peak wavelength of $600\text{ nm} \pm 10\text{ nm}$ can be used.

NOTE Addition of the buffer solution, with an initial pH of 3,8 (5.2.1), to the sample flow containing 0,025 mol/l NaOH (ratio 1,0 sample : 0,42 buffer) leads to a final pH of 4,0. Varying the pH of the distillation between 3,8 and 5,7 leads to very small fluctuations in the recovery of hydrogen cyanide (only a few percent).

4.3 Determination of easily-liberatable cyanide content

Before distillation, and without UV decomposition, a zinc sulfate solution is added to the sample flow, such that any complex iron cyanide(s) present are precipitated as zinc-cyanoferrate complexes. The hydrogen cyanide present during these conditions is separated by distillation at a temperature of 125 °C and photometrically determined as described under 4.2.

NOTE The added zinc sulfate replaces a water reagent (see Table 1) used in the determination of total cyanide content.

5 Reagents

Only use reagents of recognized analytical grade and demineralized water or distilled water according to ISO 3696, Grade 1 or 2.

5.1 General reagents

5.1.1 Hydrochloric acid $w(\text{HCl}) = 37\%$.

5.1.2 Hydrochloric acid $c(\text{HCl}) = 1\text{ mol/l}$.

Dilute 83 ml of hydrochloric acid (5.1.1) with water to 1 000 ml.

5.1.3 Hydrochloric acid $c(\text{HCl}) = 0,1\text{ mol/l}$.

Dilute 100 ml of 1 mol/l hydrochloric acid (5.1.2) with water to 1 000 ml.

5.1.4 Sodium hydroxide solution $c(\text{NaOH}) = 2,5 \text{ mol/l}$.

Dissolve 100 g of NaOH in water and dilute with water to 1 000 ml. Store in a polyethylene bottle.

5.1.5 Sodium hydroxide solution $c(\text{NaOH}) = 1 \text{ mol/l}$.

Dissolve 40 g of NaOH in water and dilute with water to 1 000 ml. Store in a polyethylene bottle.

5.1.6 Sodium hydroxide solution $c(\text{NaOH}) = 0,025 \text{ mol/l}$.

Dilute 25 ml of 1 mol/l sodium hydroxide ([5.1.5](#)) with water to 1 000 ml.

5.1.7 Detergent solution, polyoxyethylenelaurylether $\text{C}_{12}\text{H}_{25}(\text{OCH}_2\text{CH}_2)_n\text{OH}$, $n \approx 23$.

Dissolve 30 g of polyoxyethylenelaurylether by adding small quantities to 100 ml water and mix thoroughly.

NOTE This solution is commercially available as Brij-35. This information is given for the convenience of users of this document and does not constitute an endorsement by ISO of this product. To dissolve the Brij-35, the temperature can be raised to 40 °C.

5.1.8 Indicator solution

Dissolve 0,02 g of p-dimethylaminobenzylindenerhodanine in 100 ml acetone. This solution can be stored for at least one week if kept in the dark and refrigerated (5 ± 3) °C.

5.1.9 Silver nitrate solution, $c(\text{AgNO}_3) = 0,01 \text{ mol/l}$

Dissolve 1,689 7 g of silver nitrate in approximately 400 ml water and dilute to 1 litre in a volumetric flask with water. Check the actual concentration of the 0,01 mol/l silver nitrate by titration with sodium chloride according to ISO 9297 on a two weekly basis. Store this solution in the dark in a brown glass bottle. Prepare a fresh solution monthly.

5.1.10 Silver nitrate solution, $c(\text{AgNO}_3) = 0,001 \text{ mol/l}$

Prepare daily from 0,01 mol/l silver nitrate ([5.1.9](#)). Add 25,00 ml of 0,01 mol/l silver nitrate to a 250 ml volumetric flask and dilute to 250 ml with water. Store this solution in the dark in a brown glass bottle.

5.1.11 20 % V/V Denatured ethanol

Mix 100 ml denatured ethanol with 400 ml water. This solution is stable for at least one year.

5.2 Reagents for determination of cyanide**5.2.1 Distillation buffer** (pH = 3,8)

Dissolve 50 g of citric acid, $\text{C}_6\text{H}_8\text{O}_7 \cdot \text{H}_2\text{O}$ in 200 ml water. Add 120 ml sodium hydroxide solution 1 mol/l ([5.1.5](#)) and if necessary adjust the pH to 3,8 with sodium hydroxide solution 1 mol/l ([5.1.5](#)), and dilute to 500 ml with water. This solution is stable for three months if it is kept in the dark and refrigerated (5 ± 3) °C.

NOTE See [Annex C](#) for details of an alternative modified buffer to overcome flow stability problems observed with some continuous flow systems caused by the high concentration of citric acid in this reagent.

5.2.2 Zinc sulfate solution

Dissolve 10 g of zinc sulfate heptahydrate, $\text{ZnSO}_4 \cdot 7 \text{H}_2\text{O}$, in 750 ml water, mix and dilute to 1 000 ml with water. This solution is stable at least 1 year.