

SVENSK STANDARD

SS-EN ISO 10555-4:2013



Fastställt/Approved: 2013-08-09
Publicerad/Published: 2013-08-19
Utgåva/Edition: 2
Språk/Language: engelska/English
ICS: 11.040.20; 11.040.25

Sterila intravaskulära katetrar för engångsbruk – Del 4: Ballongdilatationskatetrar (ISO10555-4:2013)

Intravascular catheters – Sterile and single-use catheters – Part 4: Balloon dilatation catheters (ISO 10555-4:2013)

This preview is downloaded from www.sis.se. Buy the entire standard via <https://www.sis.se/std-98663>

Standarder får världen att fungera

SIS (Swedish Standards Institute) är en fristående ideell förening med medlemmar från både privat och offentlig sektor. Vi är en del av det europeiska och globala nätverk som utarbetar internationella standarder. Standarder är dokumenterad kunskap utvecklad av framstående aktörer inom industri, näringsliv och samhälle och befrämjar handel över gränser, bidrar till att processer och produkter blir säkrare samt effektiviserar din verksamhet.

Delta och påverka

Som medlem i SIS har du möjlighet att påverka framtida standarder inom ditt område på nationell, europeisk och global nivå. Du får samtidigt tillgång till tidig information om utvecklingen inom din bransch.

Ta del av det färdiga arbetet

Vi erbjuder våra kunder allt som rör standarder och deras tillämpning. Hos oss kan du köpa alla publikationer du behöver – allt från enskilda standarder, tekniska rapporter och standardpaket till handböcker och onlinetjänster. Genom vår webbtjänst e-nav får du tillgång till ett lättnavigerat bibliotek där alla standarder som är aktuella för ditt företag finns tillgängliga. Standarder och handböcker är källor till kunskap. Vi säljer dem.

Utveckla din kompetens och lyckas bättre i ditt arbete

Hos SIS kan du gå öppna eller företagsinterna utbildningar kring innehåll och tillämpning av standarder. Genom vår närhet till den internationella utvecklingen och ISO får du rätt kunskap i rätt tid, direkt från källan. Med vår kunskap om standarders möjligheter hjälper vi våra kunder att skapa verklig nytta och lönsamhet i sina verksamheter.

Vill du veta mer om SIS eller hur standarder kan effektivisera din verksamhet är du välkommen in på www.sis.se eller ta kontakt med oss på tel 08-555 523 00.



Standards make the world go round

SIS (Swedish Standards Institute) is an independent non-profit organisation with members from both the private and public sectors. We are part of the European and global network that draws up international standards. Standards consist of documented knowledge developed by prominent actors within the industry, business world and society. They promote cross-border trade, they help to make processes and products safer and they streamline your organisation.

Take part and have influence

As a member of SIS you will have the possibility to participate in standardization activities on national, European and global level. The membership in SIS will give you the opportunity to influence future standards and gain access to early stage information about developments within your field.

Get to know the finished work

We offer our customers everything in connection with standards and their application. You can purchase all the publications you need from us - everything from individual standards, technical reports and standard packages through to manuals and online services. Our web service e-nav gives you access to an easy-to-navigate library where all standards that are relevant to your company are available. Standards and manuals are sources of knowledge. We sell them.

Increase understanding and improve perception

With SIS you can undergo either shared or in-house training in the content and application of standards. Thanks to our proximity to international development and ISO you receive the right knowledge at the right time, direct from the source. With our knowledge about the potential of standards, we assist our customers in creating tangible benefit and profitability in their organisations.

If you want to know more about SIS, or how standards can streamline your organisation, please visit www.sis.se or contact us on phone +46 (0)8-555 523 00



Europastandarden EN ISO 10555-4:2013 gäller som svensk standard. Detta dokument innehåller den officiella engelska versionen av EN ISO 10555-4:2013.

Denna standard ersätter SS-EN ISO 10555-4, utgåva 1.

The European Standard EN ISO 10555-4:2013 has the status of a Swedish Standard. This document contains the official version of EN ISO 10555-4:2013.

This standard supersedes the Swedish Standard SS-EN ISO 10555-4, edition 1.

© Copyright/Upphovsrätten till denna produkt tillhör SIS, Swedish Standards Institute, Stockholm, Sverige. Användningen av denna produkt regleras av slutanvändarlicensen som återfinns i denna produkt, se standardens sista sidor.

© Copyright SIS, Swedish Standards Institute, Stockholm, Sweden. All rights reserved. The use of this product is governed by the end-user licence for this product. You will find the licence in the end of this document.

Upplysningar om sakinnehållet i standarden lämnas av SIS, Swedish Standards Institute, telefon 08-555 520 00. Standarder kan beställas hos SIS Förlag AB som även lämnar allmänna upplysningar om svensk och utländsk standard.

Information about the content of the standard is available from the Swedish Standards Institute (SIS), telephone +46 8 555 520 00. Standards may be ordered from SIS Förlag AB, who can also provide general information about Swedish and foreign standards.

Denna standard är framtagen av kommittén för Förbrukningsmaterial inom sjukvården, SIS/TK 330.

Har du synpunkter på innehållet i den här standarden, vill du delta i ett kommande revideringsarbete eller vara med och ta fram andra standarder inom området? Gå in på www.sis.se - där hittar du mer information.

EUROPEAN STANDARD

EN ISO 10555-4

NORME EUROPÉENNE

EUROPÄISCHE NORM

July 2013

ICS 11.040.25

Supersedes EN ISO 10555-4:1997

English Version

Intravascular catheters - Sterile and single-use catheters - Part 4: Balloon dilatation catheters (ISO 10555-4:2013)

Cathéters intravasculaires - Cathéters stériles et non réutilisables - Partie 4: Cathéters de dilatation à ballonnets (ISO 10555-4:2013)

Intravaskuläre Katheter - Sterile Katheter zur einmaligen Verwendung - Teil 4: Ballondilatationskatheter (ISO 10555-4:2013)

This European Standard was approved by CEN on 29 May 2013.

CEN members are bound to comply with the CEN/CENELEC Internal Regulations which stipulate the conditions for giving this European Standard the status of a national standard without any alteration. Up-to-date lists and bibliographical references concerning such national standards may be obtained on application to the CEN-CENELEC Management Centre or to any CEN member.

This European Standard exists in three official versions (English, French, German). A version in any other language made by translation under the responsibility of a CEN member into its own language and notified to the CEN-CENELEC Management Centre has the same status as the official versions.

CEN members are the national standards bodies of Austria, Belgium, Bulgaria, Croatia, Cyprus, Czech Republic, Denmark, Estonia, Finland, Former Yugoslav Republic of Macedonia, France, Germany, Greece, Hungary, Iceland, Ireland, Italy, Latvia, Lithuania, Luxembourg, Malta, Netherlands, Norway, Poland, Portugal, Romania, Slovakia, Slovenia, Spain, Sweden, Switzerland, Turkey and United Kingdom.



EUROPEAN COMMITTEE FOR STANDARDIZATION
COMITÉ EUROPÉEN DE NORMALISATION
EUROPÄISCHES KOMITEE FÜR NORMUNG

Management Centre: Avenue Marnix 17, B-1000 Brussels

Contents		Page
Foreword		iv
1 Scope		1
2 Normative references		1
3 Terms and definitions		1
4 Requirements		1
4.1 General.....		1
4.2 Radio-detectability.....		1
4.3 Designation of nominal size.....		1
4.4 Physical requirements.....		2
4.5 Information to be supplied by the manufacturer.....		2
Annex A (normative) Test for balloon rated burst pressure (RBP)		3
Annex B (normative) Balloon fatigue test for freedom from leakage and damage on inflation		5
Annex C (normative) Test for balloon deflation time		7
Annex D (normative) Test for balloon diameter to inflation pressure		9
Annex E (informative) Guidance on the selection of balloon materials		11
Annex ZA (informative) Relationship between this European Standard and the Essential Requirements of EU Directive 93/42/EEC amended by Directive 2007/47/EEC		12
Bibliography		1

Foreword

This document (EN ISO 10555-4:2013) has been prepared by Technical Committee ISO/TC 84 “Devices for administration of medicinal products and intravascular catheters” in collaboration with Technical Committee CEN/TC 205 “Non-active medical devices” the secretariat of which is held by DIN.

This European Standard shall be given the status of a national standard, either by publication of an identical text or by endorsement, at the latest by January 2014, and conflicting national standards shall be withdrawn at the latest by January 2014.

Attention is drawn to the possibility that some of the elements of this document may be the subject of patent rights. CEN [and/or CENELEC] shall not be held responsible for identifying any or all such patent rights.

This document supersedes EN ISO 10555-4:1997.

This document has been prepared under a mandate given to CEN by the European Commission and the European Free Trade Association, and supports essential requirements of EU Directive(s).

For relationship with EU Directive(s), see informative Annex ZA, which is an integral part of this document.

According to the CEN-CENELEC Internal Regulations, the national standards organizations of the following countries are bound to implement this European Standard: Austria, Belgium, Bulgaria, Croatia, Cyprus, Czech Republic, Denmark, Estonia, Finland, Former Yugoslav Republic of Macedonia, France, Germany, Greece, Hungary, Iceland, Ireland, Italy, Latvia, Lithuania, Luxembourg, Malta, Netherlands, Norway, Poland, Portugal, Romania, Slovakia, Slovenia, Spain, Sweden, Switzerland, Turkey and the United Kingdom.

Endorsement notice

The text of ISO 10555-4:2013 has been approved by CEN as EN ISO 10555-4:2013 without any modification.

Intravascular catheters — Sterile and single-use catheters —

Part 4: Balloon dilatation catheters

1 Scope

This part of ISO 10555 specifies requirements for balloon dilatation catheters supplied in the sterile condition, and intended for single use.

2 Normative references

The following documents, in whole or in part, are normatively referenced in this document and are indispensable for its application. For dated references, only the edition cited applies. For undated references, the latest edition of the referenced document (including any amendments) applies.

ISO 594-1, *Conical fittings with a 6 % (Luer) taper for syringes, needles and certain other medical equipment — Part 1: General requirements*¹⁾

ISO 594-2, *Conical fittings with a 6 % (Luer) taper for syringes, needles and certain other medical equipment — Part 2: Lock fittings*¹⁾

ISO 10555-1, *Intravascular catheters — Sterile and single-use catheters — Part 1: General requirements*

3 Terms and definitions

For the purposes of this document, the terms and definitions given in ISO 10555-1 and the following apply.

3.1

balloon dilatation catheter

intravascular catheter fitted with a balloon near the distal end, which is introduced into an artery or vein to dilate a part or parts of the vascular system

4 Requirements

4.1 General

Unless otherwise specified in this part of ISO 10555, balloon dilatation catheters shall comply with ISO 10555-1.

4.2 Radio-detectability

The position of the balloon shall be radio detectable when the catheter has been inserted into the body.

4.3 Designation of nominal size

The nominal size of the catheter shall be designated by the following:

- a) diameter(s) expressed in millimetres of the inflated balloon(s) or, for multidiameter balloon(s), the diameter of each portion at recommended pressure;

1) Upon its publication, ISO 80369-7 will replace ISO 594-1 and ISO 594-2.

- b) effective length of the balloon at recommended pressure;
- c) diameter of the largest guidewire that can be used with the catheter, if applicable.

NOTE Where a balloon dilation catheter (see Figure B.1) is used as a stent delivery system, refer to the appropriate standard for stents for designation of nominal size.

4.4 Physical requirements

4.4.1 Balloon rated burst pressure (RBP)

Determine the burst pressure with an appropriate safety margin when tested in accordance with [Annex A](#). Longitudinal failure is the desirable balloon failure mode.

4.4.2 Balloon fatigue; freedom from leakage and damage on inflation

Evaluate the ability of the balloon to withstand repeated inflation cycles to the RBP. When tested as described in [Annex B](#), there shall be no leakage or evidence of damage, such as herniation or bursting of the catheter.

4.4.3 Balloon deflation time

Determine the time required to deflate the balloon from the RBP as described in [Annex C](#).

4.4.4 Balloon diameter to inflation pressure

Determine the relationship between the balloon diameter and the balloon inflation pressure as described in [Annex D](#).

4.5 Information to be supplied by the manufacturer

Information supplied by the manufacturer shall comply with ISO 10555-1 and shall also include the following:

- a) nominal size of the catheter, as designated in [4.3](#);
- b) position(s) of radio-detectable marker(s);
- c) RBP of the balloon, expressed in kilopascals;
- d) balloon inflation pressure, expressed in kilopascals, required to achieve the nominal balloon diameter(s);
- e) guidewire, guide catheter or sheath or introducer compatibility and size recommendations appropriate to the intended clinical use.

NOTE Units of measurement systems other than those specified in this part of ISO 10555 can additionally be used.

Annex A (normative)

Test for balloon rated burst pressure (RBP)

A.1 Principle

The purpose of this test is to determine the RBP of the balloon.

A.2 Apparatus

A.2.1 Recommended guidewire or equivalent.

A.2.2 Water bath, controlled at (37 ± 2) °C.

A.2.3 Leak detection mechanism, e.g. dye in test fluid, pressure drop monitor, flow rate monitor.

A.2.4 Fluid for inflation, e.g. room temperature water or other justified clinically relevant media.

A.2.5 Timing mechanism, with specified accuracy.

A.2.6 Pressure generating device, fitted with a means of measuring pressure with an accuracy of ± 5 % of the reported value and maintaining the inflation pressure and fitted with a male 6 % (Luer) taper, complying with ISO 594-1 or ISO 594-2 as applicable, for connection to the catheter.

A.3 Test procedure

A.3.1 Fill the pressure generating device (A.2.6) with fluid for inflation.

A.3.2 If the instructions for use specify that a guidewire should be used during balloon inflation, insert the appropriate guidewire (A.2.1) in the device.

A.3.3 Connect the pressure generating device to the catheter under test and immerse at least the whole of the balloon portion(s) in the water bath (A.2.2) at (37 ± 2) °C.

A.3.4 Allow the catheter to equilibrate for a minimum of 2 min.

A.3.5 Inflate the balloon using a pre-determined pressure profile versus time until the catheter bursts or fails. Record the burst pressure, failure mode and location of the failure.

A.4 Test report

The test report shall include the following information:

- a) identity of the catheter;
- b) mean burst pressure, RBP and maximum, minimum and standard deviation of the burst data, expressed in kilopascals;